

ACHILLES TENDONITIS

AVERAGE DAYS WEAR

APPLICATION OVERVIEW

The Achilles tendon is the largest tendon in the body. It connects the calf muscles to the heel and is active during almost all activities including walking, jumping, and swimming. This dense tendon can withstand large forces, but can become inflamed and painful during periods of overuse. Pain results from inflammation (tendonitis) from overuse or trauma or a degenerating tendon (tendinosis). KT Tape helps treat this condition by increasing circulation and support.

WHAT YOU NEED

2 strips of KT TAPE 1 full 10" strip 1 full 10" strip cut in half

BEFORE YOU START

APPLY BEFORE ACTIVITY

Apply one hour before beginning activity

CLEAN SKIN

Clean dirt, oils and lotions from area

ACTIVATE ADHESIVE

After application rub tape vigorously to activate adhesive





STRIP ONE



0% STRETCH

ANCHOR: a full strip under the bottom of the arch with no stretch.



50% STRETCH

APPLY: the tape up the Achilles with 50% stretch.

FINISH: Lay the last two inches of the tape down without stretch.

STRIP TWO



80% STRETCH

ANCHOR: Anchor the middle of a half strip over the point of pain.

FINISH: Lay the last two inches of the tape down without stretch.

STRIP THREE



ANCHOR: If there is pain higher up on tendon, anchor middle of a second half strip to another point of pain.

FINISH: Lay the last two inches of the tape down without stretch.

